
Prion protein related peptide library

HiPep laboratories offer various prion-protein related peptides which had been designed and prepared. These peptides could discriminate structural changes of prion proteins and may inhibit or promote PrP^{Sc} generation. The interaction of these designed peptides with infectious prion obtained from mouse brain and infected mice brain homogenate has been investigated. Furthermore, by using our novel assay method to monitor structural conversion of recombinant bovine prion protein, we have found the responsible peptides for the structural conversion from bovine brain extracts. These peptide libraries are useful tool for research on protein-structural disease.

Prion protein [Homo sapiens]

1 manlgcwmlv lfvatwsdlg lckkrpkpgg wntggsrypg qgspggnryp pqggggwgqpp
61 hggwgqphg ggwgqphggg wgqphggwg qgggthsqwn kpskpktnmk hmagaaaaga
121 vvglggyvl gsamsrpiih fgdsyedry renmhrypnq vyyrpmdeys nqnnfvhdvc
181 nitikqhtvt tttkgenfte tdkmmervv eqmctqyer esqayyqrgs smvlfssppv
241 illisflif ivg

Reference

Kasai, K., Hirata, A., Ohyama, T., Nokihara, K., Yokoyama, T., Mohri, S. Novel assay with fluorescence-labelled PrP peptides for differentiating L-type atypical and classical BSEs, and scrapie, *FEBS Lett.* **2012**, 586, 325-329.

Characteristic differences of prions may account for the conformational diversity of the pathogenic isoform of prion protein (PrP^{Sc}). We applied a protein detection procedure by using fluorescent labeled peptides for detecting PrP^{Sc}. Five prion protein (PrP) related peptides were found to change significantly their fluorescent intensities with prion-affected animal samples. Their reactivity was different among atypical L-BSE, classical BSE and scrapie. The pull-down assay revealed that they precipitated PrP^{Sc} specifically. These findings suggest that fluorescent intensity changes depend on peptide PrP^{Sc} binding. This novel approach may distinguish the fine structural differences in PrP^{Sc}, which were not detected by the pull-down assay.

Hirata, A., Yajima, S., Yasuhara, T., Nokihara, K. Structural Conversion Rate Changes of Recombinant Bovine Prion by Designed Synthetic Peptides, *Int. J. Pept. Res. Ther.* **2012**, 18, 1-9.

An understanding of structural changes and self-assembly of proteins are thought to involve specific peptide-peptide interactions and it will contribute to the development of therapeutic agents and diagnosis for the detection of conformational diseases such as bovine spongiform encephalopathy or human Creutzfeldt-Jakob disease. We hypothesize that certain peptides may contribute to the conformational change of prion proteins. In contrast to the known effects of beta-sheet breaker peptides, we suggest a possible reverse phenomenon may exist in which certain peptides may induce or enhance putative beta-sheet structures. Hence designed human prion protein related peptides may be useful tools for the investigation of this hypothesis. Our paper describes the discovery of prion-related peptides which influence structural conversion of recombinant bovine prion protein.

Nokihara, K., Yajima, S., Hirata, A., Sogon, T., Yasuhara, T. Characterization of peptides obtained from digests of bovine brain which accelerate structural conversions of the recombinant bovine prion protein, *FEBS Lett.* **2013**, 587, 673-576.

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
Hirata, A., Yajima, S., Yasuhara, T., Nokihara, K. Structural Conversion Rate Changes of Recombinant Bovine Prion by Designed Synthetic Peptides *Int. J. Pept. Res. Ther.* **2012**, 18, 1-9.

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Nokihara, K., Yajima, S., Hirata, A., Sogon, T., Yasuhara, T. Characterization of peptides obtained from digests of bovine brain which accelerate structural conversions of the recombinant bovine prion protein, *FEBS Lett.* **2013**, 587, 673-676.

Nokihara, K., Yajima, S., Hirata, A. Design and syntheses of peptides which induce or enhance structural changes of recombinant bovine prion protein (rbPrP) and discovery of peptides from bovine brain which accelerate structural conversions of rbPrP, *Prion.* **2014**, 8, 117-118.

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50 ￥50,000/peptide~

all ca ￥30,000/peptide~

Each peptide is 1 mM in DMSO 100 mL stock solution

*1.5 mL microtube

Upon request tailor made sequence; Lyophilized peptides; aliquoted in a 96 well-plate can be available.

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*We have developed novel purification methods for structured peptides, such as peptides having difficult sequences. Refer also Column Information www.hipep.jp/eng.

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